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Density Standard Liquids Instructions

※Note

- Use our Standard Liquids within the following periods from date of calibration.
- Immediately use up Standard Liquids once opened.

Validity	Standard Liquids
12 months	Pure water/Iso-octane/Dichlorotoluene /Visco Liquid 5 or 50
6 months	Bromobenzene/Visco Liquid 20 or 100 or 200 or 500 or 2000/Ethanol
1 month	Sucrose Liquids

【 Cautions in handling 】

⚠ WARNING!

When handling the density standard liquids other than pure water, wear safety glasses, gloves or protective mask if necessary, and well ventilate the room.

Danger of injury on your skin or in the eyes by splashing chemicals. Also, the operator's windpipe may get hurt if toxic gas is breathe in.

⚠ WARNING!

Care should be taken on the handling of the pointed part of the aluminum seal or the tip of the needle.

The operator may cut or sting himself/herself.

⚠ WARNING!

Some density standard liquids are applicable for combustible materials, so keep from heat or flame.

The operator may burn himself/herself.

⚠ WARNING!

Dispose the density standard liquids according to the local disposal regulations or commission the disposal to an authorized industrial-waste disposal subcontractor.

The density standard liquids may have the adverse effect on the environment.


⚠ CAUTION!

Do not use the Density Standard Liquids for other than intended lab use.


Density Standard Liquids, Viscosity Liquids

- 1) Rolling up the arrow part of the aluminum sealing with a fine screwdriver or the tip of a nipper, remove the aluminum sealing and the rubber packing under it.


- 2) Collect the Liquids into a syringe

※Viscosity Liquids: Directly suck in with a syringe.


※Other Liquids: Put an aspirating needle to a syringe to suck in. Once completed, leave the needle in the vial and pull out the syringe.

- 3) Wash the syringe with the standard liquid

※Cover the tip of the syringe with paper towel to avoid the liquid splash and reciprocate the cylinder of the syringe for 20 times to clean the inside of the syringe. Discard the standard liquid used for cleaning the syringe into a waste pot.

- 4) Suck in a liquid again in the same manner as above 2). Then inject the liquid into the instrument you wish to calibrate. Do not allow entry of air bubbles while injecting.

Ethanol, Sucrose

- 1) Unpack the standard liquid

※Life the arrow marked part of the aluminum seal with a screw driver or a nipper and peel it up, then get it off and take internal gum away. Cut the inner silicone packing in a cross shape with a small blade or a razor.

※Remove the cap and the inside plug. Then put the sampling tube into the bottle and start measurement. If the sampling tube is not clean enough, the Density Standard Liquid may be contaminated. This could result in an error in measurement.

- 2) Collect the standard liquid.

※Put the clean sampling tube into the cross-section.

※ Start measurement. Ensure that no air bubbles enter in the measuring cell when injecting 8mL of the Density Standard Liquids.